

**MEDIPLUGS OPERATING ROOM
SHOES, NON-PERFORATED, 35, NAVY
BLUE**

medi Plogs


Ref. WSAM_40V

CE

 17/09/2024 13:36:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. WSAM_40V

=====

Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

Features

Bar code	3661040100036	HS Code	64029993
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	22
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	11
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	31
Packed unit / Weight	1.16	Bare Unit / Weight (kg)	1
SPECIFIC STORAGE CONDITIONS			
OTHER FEATURES			
Colors	Green	Packed unit / Quantity	1
Shoe Size	40		

17/09/2024 13:36:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr