

FOERSTER FORCEPS STRAIGHT 24 CM

Ref. IP43424



25/11/2024 20:37:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

FOERSTER FORCEPS STRAIGHT 24 CM



Ref. IP43424

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040025421	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	24
Packed unit / Weight	0.078	Bare Unit / Weight (kg)	0.077
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/11/2024 20:37:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr