

MURPHY-PEAN FORCEPS STRAIGHT 14 CM

Ref. IP21414



25/08/2024 20:19:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MURPHY-PEAN FORCEPS STRAIGHT 14 CM



Ref. IP21414

Used during surgical procedures (stop bleeding, hold compress, etc.)

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

Features

| | | | |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040023625 | HS Code | 90189090 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 23 | Bare Unit / Width (cm) | 8.5 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 0.4 |
| Packed unit / Length (cm) | 9 | Bare Unit / Length (cm) | 14 |
| Packed unit / Weight | 0.04 | Bare Unit / Weight (kg) | 0.04 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

📅 25/08/2024 20:19:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr