

## DISSECTING SET

Ref. IT01600



25/08/2024 20:33:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DISSECTING SET



Ref. IT01600

### Complete dissection kit

### Description

Consisting of:

- IP26410 Splinter forceps, 10 cm
- IP26014 Dissection forceps without teeth, 14 cm
- IS02414, Fluted probe, 14 cm
- ID07601, Dissecting needle with handle, 14 cm
- ID08601, Dissecting cutting needle with handle, 14 cm
- IC11414 Straight round-ended scissors, 14 cm
- IC09411 Straight iridectomy scissors, 11.5 cm
- ID05404, no.4 scalpel handle

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040031569	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	25	Bare Unit / Width (cm)	13.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	15	Bare Unit / Length (cm)	18
Packed unit / Weight	0.282	Bare Unit / Weight (kg)	0.282
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			

Packed unit / Quantity 45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr