


STETHOSCOPE MAGISTER II DUAL, SILVER



Ref. 507252



CE



28/11/2024 23:13:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

STETHOSCOPE MAGISTER II DUAL, SILVER



Ref. 507252

=====

Description

=====

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446042134	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Packed unit / Weight	0.43
Bare Unit / Weight (kg)	0.105		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Coral	Packed unit / Quantity	1

28/11/2024 23:13:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr