

AUDIOMETER SIBELMED SIBELSOUND DUO A OPTION AUTO HUGHSON WESTLAKE



Ref. DSI009133



17/09/2024 18:54:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

AUDIOMETER SIBELMED SIBELSOUND DUO A OPTION AUTO HUGHSON WESTLAKE



Ref. **DSI009133**

A new concept of audiometry for occupational medicine, general care and schools.

Description

Portable and lightweight, the SIBELSOUND DUO audiometer is the ideal device for performing audiometry tests in any environment.

Performs audiometry tests in automatic mode, which requires minimal technician intervention.

Follows the Hughson-Westlake automatic audiometry protocol.

Optional: export test files in standard formats (PDF, XML and CSV).

Features

CE Marking	CE0318	Medical device	Class IIa
Bar code	3661041010754	HS Code	90181990
Warranty (months)	12	Country	Spain
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	32	Packed unit / Height (cm)	14
Packed unit / Length (cm)	44	Packed unit / Weight	2.7
Bare Unit / Weight (kg)	2.7		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

17/09/2024 18:54:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr