

FLEXIBLE CHILD DIGITAL SPO2 SENSOR (90CM CABLE)

SUNTECH

Ref. DSU250904

CE
0413



 30/08/2024 09:34:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

FLEXIBLE CHILD DIGITAL SPO2 SENSOR (90CM CABLE)

SUNTECH

Ref. **DSU250904**

Description

Flexible child digital SpO2 sensor (90 cm cable)

Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	8698870322296	HS Code	90181990
Country	Turkey		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	2	Packed unit / Height (cm)	16
Packed unit / Length (cm)	14	Packed unit / Weight	0.01
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 09:34:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr