

CHROME-PLATED LAUBRY CLINIC / YELLOW



Ref. 507349



CE



30/08/2024 18:25:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CHROME-PLATED LAUBRY CLINIC / YELLOW



Ref. 507349

The famous cardiophone, vintage look for a stethoscope of very high precision

Description

Stethoscope suitable for liquid decontamination. Equipped with an ear tip closure system to ensure water tightness. It is used in isolation kits and helps prevent cross-contamination.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009182	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.415
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:25:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr