

NON CHILL RING (BLACK) FOR DUAL PULSE II SPENGLER




Ref. 507946B



 30/08/2024 18:51:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

NON CHILL RING (BLACK) FOR DUAL PULSE II SPENGLER



Ref. 507946B

Description

Spare part for SPENGLER DUAL PULSE® II stethoscopes.
Black anti-cold ring for large roof.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030667	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.002
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:51:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr