

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER



Ref. 221051



02/12/2024 19:53:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER



Ref. 221051

For SpO2 spot check.

Description

For SpO2 spot checks.
 Easy to use. Simply press the button to switch on.
 Reading comfort measures thanks to the rotatable TFT 2 colors screen :

- Simultaneous display in less than 10 sec. :
- Pulsed oxygen saturation (% SpO2),
- Heart rate,
- The pulse strength (bar graph),
- Plethysmography curve.- Low battery indicator.

- Customizable Display: 6 display modes.
 - 5 brightness levels. Switches off automatically after 5 seconds when not in use.
 Comes with wrist-strap, cover and 2 batteries.

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446025465	HS Code	90181910
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	4
Packed unit / Length (cm)	11	Packed unit / Weight	0.13
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/12/2024 19:53:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr