

Ref. E1.417.012



03/12/2024 23:29:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

Ref. E1.417.012

LUXAMED® laryngoscope, the quality for all applications!

Description

Laryngoscopes are used for intubation in operating room and in every emergency cases such as respiratory emergency and resuscitation. They are also used to perform laryngoscopies for therapeutic purpose.

The main part of the LUXAMED® laryngoscope is the handle. Its quality is especially characterized by illumination and longevity.

Lateral grooves ensure easy cleaning and secure grip in the hand.

They are available with different types of blades: Macintosh, Miller, and McCoy with different sizes.

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960348	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	45
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

 03/12/2024 23:29:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr