

LUXASCOPE OPHTHALMOSCOPE LED 3.7 V - LUXAMED

Ref. F1.426.914




Rechargeable USB



11/12/2024 21:16:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LUXASCOPE OPHTHALMOSCOPE LED 3.7 V - LUXAMED

Ref. F1.426.914

Direct ophthalmoscope in pocket size

Description

Direct ophthalmoscope in pocket size

Correction values: + 20 bis - 20 D (± 1 2 3 4 6 8 10 15 20 D)

5 different apertures: reticle, large circle, small circle, semicircle, green filter

Illumination: LED illumination

Can be dimmed for an individual glare-free brightness

Automatic switchoff after 3 minutes

Life expectancy of LEDs: approx. 50.000 hours by constant current drive

Colour temperature: approx. 4.200 K

Optimized design of head

Two-component handle (aluminium/fibre-glass reinforced), also combinable with LuxaScope LED otoscope- and dermatoscope head

Electronic ON/OFF push-button switch

Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive

Incl. recharger with USB plug (adapter EU/UK/US)

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964346	HS Code	90185090
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	White		

11/12/2024 21:16:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr