

## SUCTION UNIT CA-MI ASKIR 36 BR BOTTLE / BAG UU FLOVAC 2L



Ref. BCASPI29



13/12/2024 22:21:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SUCTION UNIT CA-MI ASKIR 36 BR BOTTLE / BAG UU FLOVAC 2L



Ref. BCASPI29

### Description

Mobile surgical Suction Unit designed to drain body fluids and perform oral, nasal and tracheal aspiration in adults and children.

- Designed for professional use and equipped with a powerful suction system that can drain up to 36 l/min max, Askir 36BR is lightweight and easy to carry.

- Three different operating options: The AC/DC adapter, rechargeable battery and 12V unit make this surgical Suction Unit a versatile unit guaranteeing CONTINUOUS interventions without the risk of overheating.

- Equipped with a visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.

- Optional roller foot and carrying case•

Supplied accessories:

- Fluid collection jar with an anti-reflux valve with different options

- 8 x 14 Silicone hoses (autoclavable) and conical connectors o 10-11-12 mm

- Antibacterial & Hydrophobic filter (single patient)

- CH20 cannula

- Power cord with Schuko plug, 12V automotive adapter, universal AC/DC adapter

- Technical characteristics:

- Motor: Oil-free piston pump requiring no maintenance

- Power supply: 14V 4A with Universal AC / DC adapter, rechargeable internal battery (12V 4A) - Automotive adapter (12V 4A)

- Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

- Max free air flow: 36 l/min

- Noise level: 63 dBA


- Weight: 4.9 kg

- Dimensions: 35 x 21 x 18 cm

 13/12/2024 22:21:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SUCTION UNIT CA-MI ASKIR 36 BR BOTTLE / BAG UU FLOVAC 2L



Ref. BCASPI29


### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	8054619100394	HS Code	90189084
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	26	Packed unit / Height (cm)	30
Packed unit / Length (cm)	40	Packed unit / Weight	5.6
Bare Unit / Weight (kg)	4.39		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 13/12/2024 22:21:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)