

ASKIR 118 WM (WALL MOUNTED) SUCTION UNIT, 1L JAR IN MAKROLON CA-MI



Ref. BCASPI306



13/12/2024 21:58:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ASKIR 118 WM (WALL MOUNTED) SUCTION UNIT, 1L JAR IN MAKROLON CA-MI



Ref. **BCASPI306**

Electrical medical device for nasal and tracheal aspiration of body fluids in children or adults orally. It has been studied specifically for use in ambulance cars or in an emerging environment.

Description

Main Applications: Car ambulance, Emergencies, Rescues

- Large LCD display for easy reading of vacuum values with vacuum adjustment keys to increase suction accuracy.
- Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / 8x14 Silicone hoses (autoclavable) and conical connector 10-11-12 mm /
- Antibacterial & Hydrophobic filter (single patient) + Cannula CH20 / Wall mount for ambulance car (Dimensions cm 20x16x16h)
- Long battery life and low noise level during operation.
- Power supply: Internal rechargeable Li-Ion battery (14.8V 5.2A) - Autonomy 70 min - charging time 360 min / Wall mount for ambulance car (12V 4A)
- PROXIMITY function is used to activate or deactivate the device without touch contacts and avoids cross-contamination between two patients.
- Long-lasting brushless motor.
- No smell of coal residues
- Maximum vacuum (adjustable): -75 kPa -563 mmHg -0.75 Bar
- Max free air flow: Hg 26 l/m (without jar)
- Operating cycle: Uninterrupted (connected to the mains)
- Weight: 2.7 kg / Dimensions: 35 x 19 x 15 cm
- Carrying case included

 13/12/2024 21:58:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ASKIR 118 WM (WALL MOUNTED) SUCTION UNIT, 1L JAR IN MAKROLON CA-MI



Ref. [BCASPI306](#)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	8054619100646	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	35
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	19
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	5.2	Bare Unit / Weight (kg)	2.7
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

13/12/2024 21:58:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr