

## ASKIR C30 BR POLYCARBONATE JARS 2X2L SUCTION UNIT



Ref. BCASPI30



13/12/2024 22:09:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

# ASKIR C30 BR POLYCARBONATE JARS 2X2L SUCTION UNIT



Ref. BCASPI30

Intended for professional use.

## Description

Ca-mi Askir C30 BR Suction Unit, autoclavable polycarbonate jars, 2 x 2000 ml

For professional use

Mobile surgical Suction Unit for draining body fluids and for oral, nasal and tracheal aspiration in adults and children.

Equipped with a powerful suction system that can drain up to 36 l/min max., two collection jars and a system on five casters, making it easy to transport.

Main applications: Emergency Department, Endoscopy, Gastroenterology, General Surgery, Obstetrics, Dentistry.

- Two different operating options: The AC/DC adapter and rechargeable batteries make this surgical Suction Unit an approved multipurpose unit, ensuring continuous interventions without the risk of overheating.

- Equipped with a visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.

-

Accessories provided: 2 fluid collection jars with an anti-reflux valve with different options / 8 x 14 Silicone hoses (autoclavable) and conical connectors 10-11-12mm / Antibacterial & hydrophobic filters (single patient) / CH20 cannula / Universal AC/DC adapter (100-240V 50/60Hz)-

Engine: Oil-free piston pump requiring no maintenance

- Power supply: 14V 4A with Universal AC / DC adapter, Rechargeable Internal Battery (12V 4A) / Autonomy 60min – Charging time 240min

-

Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 36 l/min

-

Noise level: 65 dBA

-

Operating cycle: uninterrupted (connected to the mains)

-

Weight: 7.89 kg-

Dimensions: 32 x 30 x 99 cm

 13/12/2024 22:09:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## ASKIR C30 BR POLYCARBONATE JARS 2X2L SUCTION UNIT



Ref. BCASPI30

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139166	HS Code	90189084
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	11.6
Bare Unit / Weight (kg)	7.89		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			

13/12/2024 22:09:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)