


SUCTION UNIT CA-MI NEW HOSPIVAC 350, FLOVAC SINGLE- USE POUCHES, 2 X 3000ML-



Ref. BCASPI51



13/12/2024 21:55:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 350, FLOVAC SINGLE- USE POUCHES, 2 X 3000ML-



Ref. **BCASPI51**

Description

-Basic version 3

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Available in 3 versions: Basic, FS or full• Accessories for the Basic version:- Fluid collection jar with an anti-reflux valve

8 x 14 Silicone hoses (autoclavable) and conical connectors o 10-11-12 mm

-
Antibacterial & hydrophobic filter (single patient)

-
Power cord with Schuko

plug• Accessories for the FS version:

-
Fluid collection jar with an anti-reflux

valve- 8 x 14 Silicone pipes (autoclavable) and conical connectors o 10-11-12 mm- Antibacterial & hydrophobic filter (single patient)

-
Power cord with Schuko plug

-
Pedal

• Accessories for the Full version:

-
Fluid collection jar with an anti-reflux valve

-
8 x 14 Silicone hoses (autoclavable) and conical connectors o 10-11-12 mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug

-
Pedal

- Electronic Transfer System between collection jars

• Technical characteristics

:- Engine: Oil-free piston pump requiring no maintenance

-
Power supply: 230V-50Hz

13/12/2024 21:55:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 350, FLOVAC SINGLE- USE POUCHES, 2 X 3000ML-



Ref. [BCASPI51](#)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139371	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	13
Bare Unit / Weight (kg)	13		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 13/12/2024 21:55:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr