

## CA-MI NEW HOSPIVAC 400 SUCTION UNIT, AUTOCLAVABLE JARS, 2 X 4000 ML-FS 4 VERSION,



Ref. BCASPI58



13/12/2024 21:57:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW HOSPIVAC 400 SUCTION UNIT, AUTOCLAVABLE JARS, 2 X 4000 ML-FS 4 VERSION,



Ref. BCASPI58

=====

### Description

with foot control

Designed to drain the patient's body fluids, gases, tissues and bones during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Available in 3 versions: Basic, FS or full• Accessories for the Basic version:- Fluid collection jar with an anti-reflux valve

8x14 Silicone hoses (autoclavable) and conical connectors o 10-11-12mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko

plug• Accessories for the FS version:

- Fluid collection jar with an anti-reflux

valve- 8x14 Silicone hoses (autoclavable) and conical connectors o 10-11-12mm- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug

- Pedal

• Accessories for the Full version:

- Fluid collection jar with an anti-reflux valve

- 8x14 Silicone hoses (autoclavable) and conical connectors o 10-11-12mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug,

- Pedal

- Electronic Transfer System between collection jars

Technical characteristics 13/12/2024 21:37:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW HOSPIVAC 400 SUCTION UNIT, AUTOCLAVABLE JARS, 2 X 4000 ML-FS 4 VERSION,



Ref. BCASPI58

### Features

|                             |               |                           |           |
|-----------------------------|---------------|---------------------------|-----------|
| CE Marking                  | CE0123        | Medical device            | Class IIa |
| Bar code                    | 3661040139449 | HS Code                   | 90189084  |
| Warranty (months)           | 24            | Country                   | Italy     |
| DIMENSIONS AND WEIGHT       |               |                           |           |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc        |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 20        |
| Bare Unit / Weight (kg)     | 20            |                           |           |
| SPECIFIC STORAGE CONDITIONS |               |                           |           |
| Minimum temperature         | -20           | Max temperature           | 45        |
| Minimum humidity            | NA            | Maximum humidity          | 85        |
| OTHER FEATURES              |               |                           |           |
| Packed unit / Quantity      | 1             |                           |           |

13/12/2024 21:57:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)