

IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073_RC



16/12/2024 21:24:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

IDEAL+ DOUBLE HEAD ADULT, BLUE



Ref. DH7073_RC

=====

Description

- Stethoscope with double metal horn adult (O: 45/30 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in single pavilion for adults

Available in Blue, Black, Red or Pink for the pediatric version (O: 35/30 mm)

Features

| | | | |
|------------------------------------|---------------|---------------------------|------------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661041000106 | HS Code | 9018908400 |
| Warranty (months) | 12 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 12 | Packed unit / Height (cm) | 4 |
| Packed unit / Length (cm) | 25 | Packed unit / Weight | 0.15 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Colors | Pink | Packed unit / Quantity | 1 |

📅 16/12/2024 21:24:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr