

Ref. 332012



CE



19/12/2024 21:06:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr



Ref. 332012

Description

Cotton cuff with highly resistant rubber bladder.
Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.
Washable at 30-40°C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028480	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.21
Bare Unit / Weight (kg)	0.21		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/12/2024 21:06:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr