

Ref. 332049



CE



19/12/2024 21:18:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)



Ref. 332049

## Description

Cotton cuff with highly resistant rubber bladder.  
Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.  
Washable at 30-40°C.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446008772	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	2
Packed unit / Length (cm)	13	Packed unit / Weight	0.034
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 19/12/2024 21:18:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)