

## DUAL PULSE STETHOSCOPE / BLACK



Ref. 507422



CE



19/12/2024 21:31:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DUAL PULSE STETHOSCOPE / BLACK



Ref. 507422

=====

### Description

Pearly tube model. For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaural with built-in spring.

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009212	HS Code	90181990
Warranty (months)	36		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	24	Packed unit / Height (cm)	3
Packed unit / Length (cm)	11	Packed unit / Weight	0.237
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Navy blue	Packed unit / Quantity	1

📅 19/12/2024 21:31:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)