

VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518012



19/12/2024 20:59:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518012

=====

Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

Features

| | | | |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking | CE0459 | Medical device | Class Im (measuring function) |
| Bar code | 3700446008611 | HS Code | 90189010 |
| Warranty (months) | 60 | Country | France |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 14 | Packed unit / Height (cm) | 8 |
| Packed unit / Length (cm) | 20 | Packed unit / Weight | 0.7 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 70 |
| Minimum humidity | 20 | Maximum humidity | 85 |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

📅 19/12/2024 20:59:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

🌐 www.gsh-med.fr