


## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518112



19/12/2024 21:00:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518112

=====

### Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

### Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446002077	HS Code	90189010
Warranty (months)	60	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.813
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 19/12/2024 21:00:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

🌐 www.gsh-med.fr