


POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS104

DDM



CE



 26/12/2024 21:55:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. **GBS104**

Those cuffs are autoclavable!

Description

Cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

Dimensions of the bladder in cm

Cuff infant	19 x 4 cm
Cuff child	26.5 x 5 cm
Cuff small adult	35.5 x 6 cm
Cuff adult	53 x 6 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203506	HS Code	90189084
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES

26/12/2024 21:55:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr