

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)



Ref. 520627



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3
3 ANS
GUARANTEE

📅 27/12/2024 22:46:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr


☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)



Ref. 520627



Description

HIGN ACCURACY

- Single tubing with connectic with screw
- "Watchmaking" mechanism

SHOCK RESISTANT

- Non deformable metal casing
- Shock proof foam manometer protector

COMFORT

- Ambidextrous thanks to the central decompression valve
- Large pump for quick inflation
- SPENGLER Velcro fastening

CUFF

- Coton or or nylon cuff with velcro fastening
- Sizes: **XXS** (10-15 cm), **XS** (15-22 cm), **S** (20-28 cm), **M** (26-34 cm), **L** (32-42 cm) ou **XL** (42-52 cm)
- Available in blue, grey and Cartoon
- Available in Clinic® to prevent nosocomial infection
- Comes with a transport pouch

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446008864	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.691
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 27/12/2024 22:46:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr