

CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET


Ref. GCCS104



 27/12/2024 23:02:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET



Ref. GCCS104

Those cuffs are autoclavable!

Description

Cuffs and lower limb cuff for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of **1 removable silicon bladder** for easy autoclaving.

Those ones are designed in a **conical form**

Several sizes are available

Dimensions of the bladder in cm

Lower limb cuff adult	69 x 9 cm
Lower limb cuff adult L	81 x 9 cm
Lower limb cuff adult XL	85,5 x 12 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203551	HS Code	90189084
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

27/12/2024 23:02:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr