

## CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET

Ref. GCCS105



 27/12/2024 22:35:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET

Ref. GCCS105

Those cuffs are autoclavable!

### Description

Cuffs and lower limb cuff for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of **1 removable silicon bladder** for easy autoclaving.

Those ones are designed in a **conical form**

Several sizes are available

Dimensions of the bladder in cm

Lower limb cuff adult	69 x 9 cm
Lower limb cuff adult L	81 x 9 cm
Lower limb cuff adult XL	85,5 x 12 cm

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203568	Warranty (months)	12
Country	France		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

 27/12/2024 22:35:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)