

## LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED®

Ref. F1.416.974



CE



22/01/2025 22:54:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

**LUXAMED®**

Ref. F1.416.974

**Direct ophthalmoscope in pocket size**

### Description

Direct ophthalmoscope in pocket size

Refraction balance: + 20 bis - 20 D (± 1 2 3 4 6 8 10 15 20 D)

5 different apertures: reticle, large circle, small circle, semicircle, green filter

Illumination: LED

Can be dimmed for an individual glare-free brightness

Automatic switchoff after 3 minutes

Life expectancy of LEDs: approx. 50.000 hours by constant current drive

Colour temperature: approx. 4.200 K

Optimized design of head

Two-component handle (aluminium/fibre-glass reinforced), also combinable with LuxaScope LED otoscope- and dermatoscope head

Electronic ON/OFF push-button switch

Power supply: 2x AAA Alkaline batteries, inclusive

### Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964421	HS Code	90185090
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Colors	fuchsia		

 22/01/2025 22:54:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)