

EXAMINATION COUCH STAINLESS STEEL (QUICK ASSEMBLING), BLUE LAVANDER



Ref. OSA003_L



14/01/2025 22:48:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

EXAMINATION COUCH STAINLESS STEEL (QUICK ASSEMBLING), BLUE LAVANDER



Ref. OSA003_L

Quick assembly system

Description

- Structure in square stainless steel tube
- Adjustable headrest with high resistance automatic rack and pinion
- Maximum load: 150 kg
- Head and seat padded with very high density foam (30 kg / m3), whose thickness is 5 cm, covered with washable plasticized fabric, fireproof (fire resistance classification: M1)
- Removable roller holder (for rolls 50 cm wide) positionable at the head or feet
- Available in black, Gray and Lavender.

Stainless steel or Epoxy: what's the difference?
 Epoxy: the economic choice. The steel is simply coated with a resin to delay corrosion.
 Stainless steel: it is steel modified to become completely STAINLESS steel, it does not rust.

And did you know?
 Today, more than 90% of stainless steel is made from recycled steel, it's good for the environment!

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040196152	HS Code	9402900000
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Bare Unit / Width (cm)	63
Packed unit / Height (cm)	68	Bare Unit / Height (cm)	70
Packed unit / Length (cm)	190	Bare Unit / Length (cm)	190
Packed unit / Weight	27.5	Bare Unit / Weight (kg)	24.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 14/01/2025 22:48:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

🌐 www.gsh-med.fr