

INFUSION STAND S/S - PLASTIC BASE -2 + 4 S/S SECURITY HOOKS



Ref. OSA0084C



10/01/2025 04:33:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

INFUSION STAND S/S - PLASTIC BASE -2 + 4 S/S SECURITY HOOKS



Ref. OSA0084C

Stainless steel rod, plastic base

Description

- Round stainless steel tube structure
- Black plastic base with 5 branches, diam. 56 cm
- 5 double casters very stable without brake
- Height adjustable from 120 cm to 215 cm
- Maximum load per hook: 2kg
- With supports 2 and 4 stainless steel hooks
- Delivered in 2 parts: tube and rod / foot with wheels

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168197	HS Code	94029000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	10	Packed unit / Height (cm)	54
Packed unit / Length (cm)	112	Packed unit / Weight	4
Bare Unit / Weight (kg)	3		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

10/01/2025 04:33:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr