

## OSCAR 2 ABP MONITOR - STANDARD

# SUNTECH

Ref. DSU230220



CE  
0413

DISPOSITIF MEDICAL  
CLASSE  
**IIa**  
MEDICAL DEVICE

GARANTIE  
**2**  
2 YEARS  
ANS  
GUARANTEE

 24/01/2025 22:39:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# OSCAR 2 ABP MONITOR - STANDARD



Ref. DSU230220

The reference for ambulatory blood pressure monitoring

## Description

- Clinically tested by the British Hypertension Society and the European Society of Hypertension, the Oscar 2 M250 offers very high precision measurements- SphygmoCor technology (optional) that provides important information on the clinical impact of arterial stiffness and pulse wave reflection, allowing an informed examination of cardiovascular risk
- Extraordinary comfort thanks to the ORBIT cuff
- Silent
- Lightweight and compact
- A mobile application so that the patient can inform his activity via Bluetooth (optional)
- Up to 48 hours of recording
- Analysis and editing software Accuwin Pro 4.0


## Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	0110840935105978	HS Code	90181990
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	19	Packed unit / Height (cm)	19
Packed unit / Length (cm)	30	Packed unit / Weight	1.815
Bare Unit / Weight (kg)	0.284		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 24/01/2025 22:39:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)