

SEMI-AUTOMATIC DEFIBRILLATOR SAVER ONE D

Saver One®

Ref. SJSAVE09



CE
0051

DISPOSITIF MEDICAL
CLASSE
III
MEDICAL DEVICE

GARANTIE
6 ANS
6 YEARS
GUARANTEE

24/01/2025 22:40:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SEMI-AUTOMATIC DEFIBRILLATOR SAVER ONE D



Ref. SJSAVE09

Semi-automatic AED defibrillator, simple and reliable for demanding operators with visual and audible operating information

Description

Saver One D semi_automatique

- Data visualization on large LCD screen in color (12 x 8 cm)
- Fast and practical, it is the ideal solution for teams of first aiders and paramedical staff
- Possibility to configure the device and view the ECG signal in monitoring mode
- Robust, small and very light (< 2 kg)
- Protection rating (IP54)
- Works with a lithium battery, Lifespan 4 years or 300 shocks- Rechargeable battery (optional)

- Daily automatic test and status indicator (bright led)
- Display screen for a text message
- Graph bar specifying the capacity of the battery
- Information button with bilingual mode
- Integrated metronome for better assistance with cardiac massage
- USB connection and Irda port for printing on thermal printer
- Internal memory 1 Gb

Each defibrillator is delivered with: 1 waterproof and shockproof carrying case, 1 pair of adult electrode pre-connected (lifespan 28 months), 1 first aid kit (*), 1 complete set of wall signs and 1 lithium battery

(*) Contents of the rescue kit: compresses, 1 razor, 1 mouth-to-mouth face shield, 1 pair of sterile gloves, 1 pair of scissors

NOTE: European Directive 93/42/EEC has been replaced by European Regulation 2017/745 since April 5, 2017. Formerly classified as IIB, AEDs are now classified as Class III. Certificates issued before May 25, 2017 remain valid until the end of the period indicated on these certificates, namely February 2023 for our Saver One AEDs.

 24/01/2025 22:40:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SEMI-AUTOMATIC DEFIBRILLATOR SAVER ONE D



Ref. **SJSAVE09**

Features

CE Marking	CE0051	Medical device	Class III
Bar code	3661040115436	HS Code	90189084
Warranty (months)	72	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	25.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	16
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	32
Packed unit / Weight	3	Bare Unit / Weight (kg)	2.3
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

24/01/2025 22:40:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr