

## SEMI-AUTOMATIC DEFIBRILLATOR - PRO-200J - SAVER ONE PRO - LCD SCREEN

Saver One®

Ref. SJSAVE10



CE  
0051

DISPOSITIF MEDICAL  
CLASSE  
**III**  
MEDICAL DEVICE

GARANTIE  
**6**  
6 YEARS  
GUARANTEE

 24/01/2025 22:39:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SEMI-AUTOMATIC DEFIBRILLATOR - PRO-200J - SAVER ONE PRO - LCD SCREEN



Ref. **SJSAVE10**

**Professional defibrillator, versatile with 2 modes: manual or automated**

### Description

Saver One P: PROTOCOL from 50 to 200 joules

Reliable for all scenarios, it switches very easily from semi-automatic mode to manual mode if necessary. In manual mode, it offers professionals the operational capability of ACLS (Advanced cardiac life support), allowing them to select the energy level and choose the moment for shock, perform the synchronized cardioversion medical procedure, in ECG monitoring mode and configure access to device settings.

Moreover, if necessary, it can act as an AED providing semi-automatic operating capability, with visual and audible aids, for independent BLS-trained first aiders.

- Data visualization on large LCD color screen (12 x 8 cm)
- Robust, small and very light < 2 kg
- Protection rating (IP54)
- Lithium battery, life 5 years or rechargeable battery (power version 360 joules)- Daily automatic test and status indicator (bright luminous led)
- Display screen for a text message
- Graph bar specifying the capacity of the battery
- Information button with bilingual mode
- Integrated metronome for better cardiac massage assistance
- USB connection and Irda port for printing on thermal printer
- Internal memory 1 Gb

Comes with: 1 waterproof and shockproof carrying case, 1 pair of adult electrode (26 months warranty), 1 first aid kit (\*), 1 complete set of wall signs and 1 lithium battery (standard version) or rechargeable battery (power version)

(\* ) Contents of the first aid kit: compresses, 1 razor, 1 mouth-to-mouth face shield, 1 pair of sterile gloves, 1 pair of scissors

NOTE: European Directive 93/42/EEC has been replaced by European Regulation 2017/745 since 5 April 2017. Formerly classified as IIB, AEDs are now classified as Class III. Certificates issued before May 25, 2017 remain valid until the end of the period indicated on these certificates, namely February 2023 for our Saver One AEDs.

 24/01/2025 22:39:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SEMI-AUTOMATIC DEFIBRILLATOR - PRO-200J - SAVER ONE PRO - LCD SCREEN



Ref. **SJSAVE10**

### Features

CE Marking	CE0051	Medical device	Class III
Bar code	3661040128726	HS Code	90189084
Warranty (months)	72	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.5
Bare Unit / Weight (kg)	2		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

24/01/2025 22:39:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)