

DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM



Ref. GBDM102



03/02/2025 21:16:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM



Ref. **GBDM102**

Monobloc straight arm cuff with one bladder and a ribbon

Description

Monobloc PU decontaminable arm cuff are available in straight design. Made of a single monobloc unit with one bladder and a ribbon to secure the Velcro attachment.

They are available in a range of different sizes easily indentifiable by color coding of the tubing. (A : bladder width, B: bladder length, C: total length)

		A	B	C
• New born arm cuff	White tube	3.5	29	38
• Child arm cuff	Pink tube	4.5	35	48
• Small adult arm cuff	Parme tube	6	46	61
• Adult arm cuff	Sky-blue tube	8	46	67
• Large adult arm cuff	Turquoise tube	10	47	70

Monobloc arm cuffs can not be sterilized. They can be dry cleaned or immersed paying attention to maintain the tube connector out of water.

 03/02/2025 21:16:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM



Ref. **GBDM102**

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468207917	Warranty (months)	24
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

 03/02/2025 21:16:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr