

## REPLACEMENT INFRARED EAR THERMOMETER FOR NC3/NC5

# COMEN


Ref. BCTHER01



 22/01/2025 23:07:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# REPLACEMENT INFRARED EAR THERMOMETER FOR NC3/NC5

# COMEN

Ref. BCTHER01

The wireless infrared thermometer that gives a measurement in 2 to 3 seconds!

## Features

CE Marking	NONE	Medical device	No
Bar code	3661040219349	HS Code	9025199010
Warranty (months)	6	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 22/01/2025 23:07:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)