

PEAK FLOWMETER

Ref. L36147



25/08/2024 21:01:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PEAK FLOWMETER



Ref. L36147

Description

Un débitmètre offert pour 5 boîtes d'embouts

Features

CE Marking	CE0197	Medical device	Class Im (measuring function)
Bar code	3700446001780	HS Code	90261081
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.075
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	10	Maximum humidity	95
OTHER FEATURES			
Code LPP	6115365	Packed unit / Quantity	1

25/08/2024 21:01:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr