

SET OF SILICONE TUBES Ø 6X10 FLOVAC + CONNECTOR ASPIRET, ASKIR 20/30



Ref. BCASPI111



25/08/2024 20:39:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SET OF SILICONE TUBES Ø 6X10 FLOVAC + CONNECTOR ASPIRET, ASKIR 20/30



Ref. BCASPI111

Description

Set of pipes containing:- 1 pipe diam 6 mm / ext 10mm x length

140cm

- 1 pipe diam 6 mm / ext 10mm x length 26 cm

- Connec connectors 8-9-10 mm

Features

CE Marking	CE0123	Medical device	No
Bar code	3661040160818	HS Code	90330090
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	18.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	11
Packed unit / Weight	0.05	Bare Unit / Weight (kg)	0.05
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 20:39:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr