

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AUTOCLAVABLE JARS, 2 X 5000 ML- FULL 5 VERSION



Ref. BCASPI506



 13/12/2024 22:04:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AUTOCLAVABLE JARS, 2 X 5000 ML- FULL 5 VERSION



Ref. **BCASPI506**

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

Description

-with foot control and automatic suction inverter

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Accessories included:- Autoclavable

jars with overflow valve

-8x14 Autoclavable Silicone Tubes + Conical connector o 10-11-12mm

-Antibacterial and hydrophobic filter

-Safety container 220 ml

- Technical characteristics

:- Engine: Oil-free piston pump requiring no maintenance

-

Power supply: 230V-50Hz

-

Maximum vacuum (adjustable): -0.90 Bar -90 kPa -675 mmHg

- Max free air flow: 60 l/min

- Noise level: 52 dBA

- Weight: 13 kg

- Dimensions: 46 x 85 x 42 cm

 13/12/2024 22:04:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AUTOCLAVABLE JARS, 2 X 5000 ML- FULL 5 VERSION



Ref. [BCASPI506](#)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661041009932	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	13
Bare Unit / Weight (kg)	13		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

13/12/2024 22:04:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr