

## CAIX ELECTRICAL EXAMINATION COUCH

Ref. OV003\_B



29/08/2024 16:49:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## CAIX ELECTRICAL EXAMINATION COUCH



Ref. OV003\_B

Caix electric table from 560 to 960 mm - Width 670mm

### Description

- Electric couch, height adjustable from 560 to 960 mm
- Built in white epoxy rectangular tubing
- Equipped with 4 adjustable black pads, in polyamide
- Sewn upholstery with a thickness of 40 mm, foam seat and backrest with a density of 28 kg/m3
- Adjustable backrest by gas cylinder
- Equipped with a roll holder
- Foot control
- Category II motor
- Various colours to choose from, please ask us.

- Width: 670 mm
- L. 2000 × l. 670 × H. 560/960 mm
- Backrest: L. 600 × l. 670 mm
- Seat: L. 1350 × l. 670 mm
- Weight of the couch: 80 kg
- Maximum load: 220 kg

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094373	HS Code	940290000
Warranty (months)	36	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	67
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	92
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	195

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr