

ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGHT AND LARG) + POUCH, ORANGE



Ref. SBRAN08_O



23/09/2024 10:22:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGHT AND LARG) + POUCH, ORANGE



Ref. **SBRAN08_O**

Folding stretcher with balance wheel system

Description

Folding stretcher with balancing system

- High resistance aluminium structure and waterproof fabric, fire resistant and anti-cracking
- This balancing stretcher folds in length and width
- Lightweight, easy to carry
- Dimensions (L x W x H): 208 X 55 x 13 cm, unfolded; 104 x 17 x 9 cm, folded
- Weight: 8 kg
- Maximum load carried: 159 kg
- Single colour: orange
- Delivered with a stretcher bag

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094632	HS Code	94029000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	20	Bare Unit / Width (cm)	54.5
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	13.5
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	208
Packed unit / Weight	5.6	Bare Unit / Weight (kg)	4.7
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

23/09/2024 10:22:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr