

## ANTIBACTERIAL & HYDROPHOBIC FILTER FOR ASKIR C30BR




Ref. BCASPI97



 29/08/2024 15:44:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## ANTIBACTERIAL & HYDROPHOBIC FILTER FOR ASKIR C30BR



Ref. BCASPI97

### Description

- Anti-bacterial and hydrophobic filter
- Diam 64mm / Connector 11mm

### Features

CE Marking	NONE	Medical device	No
Bar code	8054619103869	HS Code	90192000
Warranty (months)	12	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	16
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	5.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	25
Packed unit / Weight	0.05	Bare Unit / Weight (kg)	0.05
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 15:44:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)