

MIXTER-O'SHAUGNESSY ARTERY FORCEPS 90° 18 CM

Ref. IP49018



09/01/2025 22:06:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

MIXTER-O'SHAUGNESSY ARTERY FORCEPS 90° 18 CM



Ref. IP49018

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040138312	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	8	Packed unit / Height (cm)	2
Packed unit / Length (cm)	19	Packed unit / Weight	0.07
Bare Unit / Weight (kg)	0.07		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

09/01/2025 22:06:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr