

ELITE BAGS - PODIA - GREY



Ref. MEPODO02B



29/08/2024 15:52:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ELITE BAGS - PODIA - GREY



Ref. MEPODO02B

New generation of this podiatry bag, with a new ergonomics more adapted to the needs of professionals.

Description

Large capacity multidisciplinary case that is perfectly adapted to the needs of professionals who need to carry a large number of instruments and devices (such as the podiatrist's micromotor)

- Possibility to organize all the equipment in one space for more efficient patient care

- Two work areas in an original main compartment, separated by a central shutter

- Detachable pocket for devices

- At the front, space for personal items

- Material: fabric 900D

- Weight: 3.20 kg

- Outdoor measures: 42 x 32 x 16 cm / 23.55 L

- Recommended max. capacity: 5 kg

- Optional trolley

- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540041	HS Code	420212
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	32
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	16
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	42
Packed unit / Weight	3.5	Bare Unit / Weight (kg)	3.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:52:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr