

ALUMINIUM ALLOY SCOOP STRETCHER



Ref. SBRAN12



29/08/2024 12:02:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr


+33(0)4 42 90 31 31

info@gsh-med.fr

ALUMINIUM ALLOY SCOOP STRETCHER



Ref. **SBRAN12**



Description

- «Spoon» function performed by the two longitudinal blades, to position the stretcher under the patient without moving him. The scoop stretcher can be used to lift, move, or transport patients.
- Length adjustable according to the patient, from 120 to 190 cm
- High resistance aluminium structure
- Dimensions (L x W x H): 120-190 x 44 x 6 cm
- Weight: 9 kg
- Maximum load carried: 159 kg

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040063669	HS Code	94029000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	45	Bare Unit / Width (cm)	44
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	9
Packed unit / Length (cm)	170	Bare Unit / Length (cm)	120
Packed unit / Weight	10	Bare Unit / Weight (kg)	8.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 12:02:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr