

RAPID CONNECTING SET FOR SPHYGMOMANOMETER HOLTEX



Ref. DS99012



 26/08/2024 19:20:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

RAPID CONNECTING SET FOR SPHYGNOMANOMETER HOLTEX



Ref. DS99012

Features

CE Marking	NONE	Medical device	No
Bar code	3661040074566	HS Code	9018199560
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	9
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14
Packed unit / Weight	0.01	Bare Unit / Weight (kg)	0.01
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 19:20:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr