

## TOBOLD METAL TONGUE DEPRESSOR WITH HANDLE 14.5 CM



Ref. IA02404



26/08/2024 19:26:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## TOBOLD METAL TONGUE DEPRESSOR WITH HANDLE 14.5 CM



Ref. IA02404

### Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016). In the "

ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

### Features

|                                    |               |                         |           |
|------------------------------------|---------------|-------------------------|-----------|
| CE Marking                         | CE            | Medical device          | Class I   |
| Bar code                           | 3661040005980 | HS Code                 | 9018.9090 |
| Warranty (months)                  | 24            | Country                 | Pakistan  |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                         |           |
| Packed unit / Width (cm)           | nc            | Bare Unit / Width (cm)  | 9         |
| Packed unit / Height (cm)          | nc            | Bare Unit / Height (cm) | 2         |
| Packed unit / Length (cm)          | nc            | Bare Unit / Length (cm) | 17        |
| Packed unit / Weight               | 0.096         | Bare Unit / Weight (kg) | 0.096     |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                         |           |
| Minimum temperature                | NA            | Max temperature         | NA        |
| Minimum humidity                   | NA            | Maximum humidity        | NA        |
| <b>OTHER FEATURES</b>              |               |                         |           |
| Packed unit / Quantity             | 1             |                         |           |

26/08/2024 19:26:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)