

HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°1/2



Ref. IS01401

Exemple de bougies Hégar double, inox.
15 tailles disponibles, nous consulter.



27/08/2024 05:46:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°1/2



Ref. IS01401

Description

Hégar Double Ended Dilator Stainless Steel

15 Available sizes, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029542	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	4.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.007	Bare Unit / Weight (kg)	0.007
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 05:46:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr