

## HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°29



Ref. IS01629



Exemple de bougies Hégar simple, inox.  
30 tailles disponibles, nous consulter



27/08/2024 06:14:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°29



Ref. IS01629

### Description

Hégar Single Ended Dilator Stainless Steel

30 Sizes available, contact us

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029993	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	3.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.118	Bare Unit / Weight (kg)	0.118
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

27/08/2024 06:14:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)