


LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

LUXAMED®

Ref. A1.616.314



Technologie CCT
2 couleurs de lumière



13/03/2025 23:19:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

LUXAMED®

Ref. A1.616.314

Innovative otoscope in pocket size

Description


- Innovative illumination: LED-Ring with 2x4 LEDs (concentric)
switchable with one button press : warm white light or cold white light,
no fibre optic,
field of vision is not restricted
- Maintenance-free, because no lamp replacement necessary
 - Automatic switchoff after 3 minutes
 - Life expectancy of LEDs: approx. 100.000 hours due to ideal operation by constant current drive
 - Colour temperature: ring 1 approx. 3.000 K, ring 2 approx. 6.000 K
 - Optimized design of head
 - Pivoting 3-fold magnifying glass with locking function
 - Insufflation port for pneumatic test
 - Two-component handle (aluminium/fibre-glass reinforced), also combinable with Luxascope Dermatoscope CCT LED head
 - Electronic ON/OFF push-button switch
 - Power supply: 2x AAA Alkaline batteries, inclusive
 - Incl. 20 disposable ear funnels, grey, each 10 x 2.5 mm and 10 x 4.0 mm

Available in :
Black : A1.616.114
Blue : A1.616.214
Grey : A1.616.314
White : A1.616.914

 13/03/2025 23:19:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

LUXAMED®

Ref. A1.616.314


Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963592	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	4
Packed unit / Length (cm)	24	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Gray		

 13/03/2025 23:19:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr