

Ref. C1.616.114



Technologie CCT
2 couleurs de lumière



13/03/2025 23:02:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

Ref. C1.616.114

Innovative dermatoscope in pocket size

Description

Innovative illumination: LED-Ring with 2x5 LEDs in concentric structure for a highly detailed image, homogenous light, switchable with one button press : warm white light or cold white light

- Removable contact plate for easy cleaning
- Maintenance-free, because no lamp replacement necessary
- Stepless zoom area with 10 x magnification for a distortion-free image
- Optimized design of head
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours due to ideal operation by constant current drive
- Intensity of illumination: 10.000 Lux
- Colour temperature: ring 1 approx. 3.000 K, ring 2 approx. 6.000 K
- Two-component-handle (aluminium/fibre-glass reinforced),
- Also combinable with LuxaScope Auris CCT LED 2.5 V otoscope-head
- Electronic ON/OFF push-button-switch
- Power supply: 2 x AAA Alkaline batteries, inclusive

Available in :

- Black : C1.616.114
- White : C1.616.914

 13/03/2025 23:02:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

Ref. C1.616.114

Features

CE Marking	NONE	Medical device	No
Bar code	04260381963615	HS Code	90189084
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		

13/03/2025 23:02:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr